HISTOSONDA c-erbB2

Referencia: CBM-0020-R5

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HISTOSONDA c-erbB2

Reference: CBM-0020-R5

Assay of one vial for 5 individual reactions. 65µl for reaction.

Product classification:

*EU countries
"For In Vitro Diagnostics"

INTRODUCTION

c-erbB2 is denominated as such as it corresponds to the human homologue of a viral gene that causes leukemia erythroblastosis in birds. This gene was isolated in a neuroblastoma/glioblastoma rat cell colony and for this reason it is also known as HER2-NEU.

This gene codes for a tyrosine kinase receptor situated in the cytoplasmic membrane and for whom the ligand is still unknown.

This protein belongs to a family of receptors whose principle member is EGF α r (epidermal growth factor receptor) also known as c-erbB1. This family is made up of four members (c-erbB1, c-erbB2, c-erbB3 and c-erbB4) that share great internal homology. Due to this homology antibodies targeted against this family of proteins can cross react with other members.

It is clinically shown that between 15 and 30% of breast cancers over express the c-erbB2 protein. This over expression implies a worse prognosis for these tumors meaning that they are more aggressive and resistant to Taxol.

As there are therapeutic monoclonal antibodies against this protein it is vital to know if a concrete tumoral case over expresses the c-erbB2 protein or not.

Until now the methods of detection have been based on immunohistochemistry with antibodies targeted against the c-erbB2 protein which generally yield confusing results and the determination of amplification of the gene that codes for this product.

Cenbimo has designed a probe for in situ hybridization that detects the transcribed RNA of the c-erbB2 gene.

The detection of messenger RNA for this protein offers multiple advantages as it eliminates the confusing aspects of immunohistochemistry, and on the other hand detects real transcription, as to see genetic amplification by FISH does not guarantee that all amplified genes are transcribing.

Histosonda c-erbB2 consists of 2 fragments of single-stranded DNA with lengths of 694 and 1143 nucleotides that are targeted against expressed sequences of this gene.

The probe incubation period is only 1 hour and as it is a chromogenic hybridisation probe the results can be viewed directly by light microscope.

INTENDED USE

For use in In Vitro Diagnosis.

Histosonda c-erbB2 is useful for the detection of c-erbB2 RNA and its associated tumors.

WARNINGS AND PRECAUTIONS

Histosonda c-erbB2 has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel.

In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

COMPONENTS

1 tube with the necessary amount of Histosonda c-erbB2.

Histosonda c-erbB2 consists of two segments of single-stranded DNA complementary to expressed RNA, with lengths of 695 and 1144 nucleotides. This probe has been labeled with Digoxigenin.

STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

INTERPRETATION OF RESULTS

Samples in which c-erbB2 expression is observed will show a brownish color in the cell cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, and conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

ASSAY LIMITATIONS

Histosonda c-erbB2 has been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

BIBLIOGRAPHY

- 1) HER2 status determination using RNA-ISH - a rapid and simple technique showing high correlation with FISH and IHC in 141 cases of breast cancer. Javier Alba¹, Javier Gutiérrez², Victoria Mary Coupe², Beatriz Fernández¹, Ángel Vázquez-Boquete¹, Jesús Alba², Jerónimo Forteza¹ and Tomás García-Caballero¹
- 2) mRNA In Situ Hybridization (Histosonda): A New Diagnostic Tool for HER2-Status in Breast Cancer. A Multicentric Spanish Study. Bernet L, Martínez Benaclocha M, Castera C, Cano Muñoz R, Sevilla F, Alba J, de Dios Barranco J, Córdoba A, García-Caballero T, Hardisson D, de Francisco Hernandez JM, Lazaro JM, Polo L, Riu F, Rezola R, Rojo F, Ruiz I, Hernández A, de la Cámara de Las Heras JM, Coupe VM.

PROCEDURE

For Digoxigenin- labeled probes

BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

HistoSonda Protocol

Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete HistoSonda Technical Manual at the email address: techserv@cenbimo.com

Before starting pre-heat the humid chamber to 62°C.

1a. Deproteinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
 - a. Xylene: 10 mins
 - b. Xylene: 5mins
 - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
 - d. Alcohol 96%: 1min X 3
 - e. Methanol containing 0.3% H₂O₂ or 3% H₂O₂ only: 5mins
3. Wash well with distilled water, leave standing in water.

1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

1. Place slides in boiling distilled water for 30 seconds and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033mg/ml in PBS pH 7,4.
2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
4. Wash well with distilled water.
5. Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment!**)

3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.

4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.!

1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS, pH 7.4 agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS, agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
8. Wash with water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in water.
10. Wash with water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: CENBIMO S.L.
C/Doctor Iglesias Otero s-n
27004-Lugo, Spain