## HISTOSONDA CITOMEGALOVIRUS (CMV)

Referencia: CBM-0018-R5

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### BIBLIOGRAFIA

- 1) The most abundantly transcribed human cytomegalovirus gene ( $\beta$ 2.7) is not essential for growth in vitro. Brian P. McSharry, 1 Peter Tomasec, 1 M. Lynne Neale 2 and Gavin W. G. Wilkinson 1
- 2) Localization of the human cytomegalovirus 2.7-Kb major early beta-gene transcripts by in situ hybridization in permissive and nonpermissive infections. TC Wu, WA Lee, MC Pizzorno, WC Au, YJ Chan, RH Hruban, GM Hutchins and GS Hayward

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## HISTOSONDA CYTOMEGALOVIRUS (CMV)

Reference: CBM-0018-R5

Assay of one vial for 5 individual reactions. 65µl for reaction.

Product classification:  
\*EU countries  
"For Research Use Only" (RUO)

### INTRODUCTION

Cytomegalovirus belongs to the β-herpes virus family and globally affects the world population. Primary infection is often a flu-like process although there are known cases of greater severity that are clinically similar to mononucleosis.

Like other herpes viruses, after primary infection CMV persists in the host permanently as a latent infection.

Immunosuppression in patients can provoke reactivation of the virus giving rise to serious lesions in a multitude of organs such as the lungs, kidneys, central nervous system, gastrointestinal apparatus, endothelium and blood cells.

Patients with ulcerative colitis can be infected by CMV and treatment guidelines modified by this situation.

After infection, CMV first begins to transcribe the so-called 'early' genes. Of these genes, one known as β2.7 (2.7 kbases) is the most abundantly transcribed and accumulates during the entire infection until it makes up 20% of all mRNA transcribed by the virus. Because of this, the gene β2.7 is an ideal candidate for CMV identification.

CENBIMO has fabricated HISTOSONDA CMV which consists of a single strand of DNA with a sequence length of 288 nucleotides. This sequence is complementary to CMV gene β2.7 mRNA. The probe DNA has been labeled with digoxigenin.

The technique of in situ hybridization presents numerous advantages over immunohistochemistry and unspecific results obtained when using antibodies can be avoided.

### INTENDED USE

For Research Use Only (RUO).

Histosonda Cytomegalovirus is intended for the detection of cells infected by CMV in any location found including: lymph node, central nervous system, retina, lung and intestine. This detection is carried out by in situ hybridization in paraffin-embedded, formalin-fixed histological sections.

### WARNINGS AND PRECAUTIONS

Histosonda Cytomegalovirus has been designed for In Vitro Diagnostic use only and must be manipulated by qualified and accordingly trained personnel.

In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

This product contains sodium azide as a preservative (NaN<sub>3</sub>). Sodium azide is below the threshold level of 0,1% (w/v) of the total, and therefore it is not considered to be a hazardous substance. Small quantities may be eliminated in the waste water system.

### COMPONENTS

1 tube the necessary amount of Histosonda Cytomegalovirus.

Histosonda Cytomegalovirus consists of a fragment of single stranded DNA of 288 nucleotides, complementary to CMV gene β2.7 mRNA. The probe DNA contains nucleotides that have been labeled with Digoxigenin.

### STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

### SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

### INTERPRETATION OF RESULTS

Samples in which Cytomegalovirus expression is observed will show a brownish color in the cell nucleus, cytoplasm or both, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample in parallel with the staining observed in positive and negative controls.

### ASSAY LIMITATIONS

Histosonda Cytomegalovirus has been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

### BIBLIOGRAPHY

- 1) The most abundantly transcribed human cytomegalovirus gene (β2.7) is no essential for growth in vitro. Brian P. McSharry, 1 Peter Tomasec, 1 M. Lynne Neale 2 and Gavin W. G. Wilkinson 1
- 2) Localization of the human cytomegalovirus 2.7-Kb major early beta-gene transcripts by in situ hybridization in permissive nonpermissive infections. TC Wu, WA Lee, MC Pizzorno, WC Au, YJ Chan, RH Hruban, GM Hutchins and GS Hayward

### PROCEDURE

For Digoxigenin- labeled probes

### BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

### HistoSonda Protocol

#### Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polimero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete **HistoSonda Technical Manual at the email address: techserv@cenbimo.com**

Before starting pre-heat the humid chamber to 62°C.

#### 1a. Deparaffinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
  - a. Xylene: 10 mins
  - b. Xylene: 5mins
  - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
  - d. Alcohol 96%: 1min X 3
  - e. Methanol containing 0.3% H<sub>2</sub>O<sub>2</sub> or 3% H<sub>2</sub>O<sub>2</sub> only: 5mins
3. Wash well with distilled water, leave standing in water.

#### 1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally **bone marrow and gastric tissues**).

1. Place slides in boiling distilled water for **30 seconds** and immediately transfer to distilled water at room temp.

**! Important:** excess heat treatment will result in background staining. Remove slides immediately. !

#### 2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033 mg/ml in PBS pH 7.4.
2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
4. Wash well with water.
5. Transfer to PBS pH 7.4 for 2mins.

**! (Important:** if the slides have been boiled only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment!**)

#### 3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.

#### 4. Washing the probe

**! The probe solution is very viscous due to large quantities of dextran. It is very important to wash the slides well to ensure the probe solution is completely removed!**

1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

#### 5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS, pH 7.4 agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS, agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
8. Wash with water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in water.
10. Wash with water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

**Manufacturer: CENBIMO S.L.**  
**C/Doctor Iglesias Otero s/n**  
**27004-Lugo, Spain**