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1. Peter J. Delves And Ivan M. Roitt: Immunoglobulin genes. In *ENCYCLOPEDIA OF IMMUNOLOGY*. Pag 1323. Second Edition, ACADEMIC PRESS LIMITED (1998).
2. Cunningham-Rundles C, Bodian C. Common variable immunodeficiency: clinical and immunological features of 248 patients. *Clin Immunol* 1999;92:34-48
3. Clark JA, Callicot PA, Brenner NA. Selective IgA deficiency in blood donors. *Am J Clin Pathol* 1983; 80:210-213.
4. Notarangelo LD, Duse M, Ugazio AG. Immunodeficiency with hyper IgM (HIM). *Immunodef Rev* 1992;3:101-121

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(DELTA OR EPSILON)

References:

Histosonda Delta Heavy Chain CBM-0007-R5
Histosonda Epsilon Heavy Chain CBM-0008-R5

Assay for 5 individual reactions in one vial. 65µl for reaction.

Product classification:
*EU countries
For In Vitro Diagnostics

INTRODUCTION

Competent B lymphocytes immunologically develop from a lymphoid stem cell until they acquire a receptor constituted of a single immunoglobulin molecule specific for a determined antigen. These immunoglobulin molecules are formed by 2 identical heavy chains and 2 identical light chains (kappa or lambda).

In humans the immunoglobulin heavy chains gene is located in the telomeric region of the long arm of chromosome 14 and spans a length of approximately 1400000 base pairs. The zones of this gene that code for the immunoglobulin heavy chain constant regions form a continuum situated at the 3' end of the gene and are ordered in the direction 5' to 3' - Cµ, Cδ, Cγ3 Cγ1, Cα1, Cγ2, Cγ4, Cε y Cα2.

In human normal or reactive lymphoid tissues the immunoglobulin producing cells appear in similar proportions to the seric distribution of immunoglobulins with the IgG cells being more abundant, then the IgA cells and finally the IgM cells, with the exception of the intestinal lamina in which the most abundant immunoglobulin is IgA. Very few plasmatic cells producing IgD are observed in lymphoid tissues with the exception of nasopharyngeal tonsils where they can be more abundant. In the B lymphocyte tumors, as a general rule, only one clone proliferates and therefore only one heavy chain is observed. The detection of monoclonality is one of the most important tools for differentiating B lymphoid tumors from reactive processes.

The "in situ" hybridization technique has an important advantage over immunohistochemistry in that there is virtually no background staining and it permits a clear visualization of the histological preparation.

INTENDED USE

For use in In Vitro Diagnosis.

The Histosondas of the Immunoglobulin Heavy Chains are useful for the study of monoclonality in lymphoid tumours, lymphoproliferative syndromes, myelomas and for the study of immunodeficiency syndromes where the chain-specific plasmatic cells are absent or diminished (supposedly one in every 350 Europeans suffer from selective deficit of IgA).

WARNINGS AND PRECAUTIONS

The Histosondas of the Immunoglobulin Heavy Chains have been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel.

In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

This product contains sodium azide as a preservative (NaN₃). Sodium azide is below the threshold level of 0,1% (w/v) of the total, and therefore it is not considered to be a hazardous substance. Small quantities may be eliminated in the waste water system.

COMPONENTS

1 tube with the necessary amount of Histosonda Delta or Epsilon for developing 5 single reactions.
The Histosondas of the Immunoglobulin Heavy Chains consist of a fragment of single-stranded DNA with a length of between 557 and 601 nucleotides, complementary to expressed RNA. These probes have been labeled with Digoxigenin.

STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study. width are sufficient to conduct the study.

INTERPRETATION OF RESULTS

Samples in which immunoglobulin heavy chain expression is observed will show a brownish color in the cell cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

ASSAY LIMITATIONS

The Histosondas of the Immunoglobulin Heavy Chains have been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Their use is not recommended for other types of samples or preparation techniques.

The correct operation of these products has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

BIBLIOGRAPHY

1. Peter J. Delves And Ivan M. Roitt: Immunoglobulin genes. In ENCYCLOPEDIA OF IMMUNOLOGY, Page 1323. Second Edition, ACADEMIC PRESS LIMITED (1998)
2. Cunningham-Rundles C, Bodian C. Common variable immunodeficiency: clinical and immunological features of 248 patients. Clin Immunol 1999;92:34-48
3. Clark JA, Callicot PA, Brenner NA. Selective IgA deficiency in blood donors. Am J Clin Pathol 1983; 80:210-213.
4. Notarangelo LD, Duse M, Ugazio AG. Immunodeficiency with hyper IgM (HIM). Immunodef Rev 1992;3:101-121

PROCEDURE

For Digoxigenin- labeled probes

BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosonda targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

HistoSonda Protocol

Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete HistoSonda Technical Manual at the email address: techserv@cenbimo.com

Before starting pre-heat the humid chamber to 62°C.

1a. Deparaffinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
 - a. Xylene: 10 mins
 - b. Xylene: 5mins
 - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
 - d. Alcohol 96%: 1min X 3
 - e. Methanol containing 0.3% H₂O₂ or 3% H₂O₂ only: 5mins
3. Wash well with distilled water, leave standing in water.

1b. Inhibition of unspecific DNA binding (optional)
Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

1. Place slides in boiling distilled water for 30 seconds and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0.033 mg/ml in PBS pH 7.4.
2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
4. Wash well with water.
5. Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment) !**

3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.

4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed!



1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS pH 7.4 agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS pH 7.4, agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
8. Wash with distilled water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in distilled water.
10. Wash with distilled water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: CENBIMO S.L.
C/Doctor Iglesias Otero s-n
27004-Lugo, Spain