HISTOSONDA DUAL KAPPA-LAMBDA

Referencia: CBM-0027-R5

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HISTOSONDA DUAL KAPPA-LAMBDA

Reference: CBM-0027-R5

Assay of one vial for 5 individual reactions. 65µl for reaction.

Product classification:
*EU countries
For In Vitro Diagnostics

INTRODUCTION

Immunologically competent B lymphocytes develop from a lymphoid stem cell until acquiring a receptor consisting of a single immunoglobulin molecule specific for a determined antigen. These immunoglobulin molecules are formed by two identical heavy chains and by two identical light chains (kappa or lambda). In humans, in normal or reactive lymphoid tissues, approximately 60% of the B lymphocytes bear kappa light chains and the remaining 40%, lambda light chains. In B lymphocyte tumors, there is generally only one clone that proliferates, be it kappa or lambda, in a monoclonal way. In humans the gene that codes for kappa light chains is located on chromosome 2, and the gene that codes for lambda light chains is located on chromosome 22. The variable regions for genetic recombination are located in the first 5' half of these genes, while the constant regions are located in the 3' half. There is only one functional gene for the constant region of kappa light chains, and four functional genes for the constant region of lambda light chains (CA1, CA2, CA3, and CA7). Detection of monoclonality is one of the most important tools for differentiating B lymphoid tumors from reactive processes.

In situ hybridization technique offers an important advantage over immunohistochemistry, as it virtually lacks background, and allows a clean and sharp viewing of the histological preparation. It is also useful to differentiate cells that have absorbed immunoglobulins, and are therefore detectable by immunohistochemistry, but in fact do not produce immunoglobulin, as occurs with the Reed-Sternberg cells of Hodgkin's disease.

INTENDED USE

For use in In Vitro Diagnosis. Histosonda Dual Kappa-Lambda is useful for the study of monoclonality in lymphoid tumors, lymphoproliferative syndromes, myelomas and for the study of immunodeficiency associated lymphoproliferative syndromes.

WARNINGS AND PRECAUTIONS

Histosonda Dual Kappa-Lambda has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel.

In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

COMPONENTS

1 tube with the necessary amount of Histosonda Dual Kappa-Lambda for developing 5 single reactions. Histosonda Dual Kappa-Lambda consists of two fragments of single-stranded DNA with lengths of 153 and 182 nucleotides that are complementary to expressed RNA. In this probe the Lambda fragment has been labeled with Digoxigenin and the Kappa fragment has been labeled with Biotin.

STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

INTERPRETATION OF RESULTS

Samples in which immunoglobulin light chains expression is observed will show color in the cell cytoplasm, which will contrast over the blue-violet background

given by hematoxylin staining. The color will be brown in the case of Lambda light chain, the color for Kappa light chains depends on the choice of chromogenic reagent used (generally red or blue). The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls for expression of this gene.

ASSAY LIMITATIONS

Histosonda Dual Kappa-Lambda has been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

PROCEDURE

Recommended protocol for Histosonda Dual probes.

BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labelled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

HistoSonda Protocol

Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP commercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete **HistoSonda Technical Manual at the email address: techserv@cenbimo.com**

Before starting pre-heat the humid chamber to 62°C.

1a. Deparaffinization

- Heat slides at 62°C for 10mins using a hot plate or incubator.
- Immerse slides in:
 - Xylene: 10 mins
 - Xylene: 5mins
 - Absolute alcohol (ethanol or isopropanol): 1min X 2
 - Alcohol 96%: 1min X 3
 - Methanol containing 0.3% H₂O₂ or 3% H₂O₂ only: 5mins
- Wash well with distilled water, leave standing in water.

1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

- Place slides in boiling distilled water for **30 seconds** and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

2. Deproteinization

- Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033mg/ml in PBS pH 7,4.
- Remove excess water from individual slides with tissue paper.
- Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
- Wash well with distilled water.
- Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5

minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment!**

3. Incubation with the probe

- Remove excess buffer from sections as described previously.
- Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
- Place the slides in a horizontal position in a humid chamber and close well.
- Incubate at 62°C for 1hr.

4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.!

- Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
- Agitate in PBS for 5mins.

5. Revealing the probe

- Dilute an Anti-Biotin (not provided with this product) to working concentration with an Anti-Digoxin antibody (not provided with this product(mouse) following the manufacturer's instructions.
- Cover slides with 100ul of the resulting primary antibody mix and incubate in a humid chamber at room temperature for 30mins.
- Wash the slides with PBS pH 7.4.
- Cover the slides (1-2drops) with secondary HRP and AP cocktail polymer (not provided with this product) attending to the type of the primary antibodies used previously and following the manufacturer's instructions. Incubate in a humid chamber at room temperature for 30mins.
- Wash slides with **TBS pH 7.4**
- Make up DAB solution (not provided with this product) following the manufacturers instructions (normally 1µl DAB in 50µl substrate for each slide)
- Apply DAB solution (50µl) and observe development by light microscope.
- Wash with **TBS pH 7.4**
- Make up Fast Red solution (not provided with this product) following the manufacturer's instructions (normally 1.5µl Fast Red in 50µl substrate for each slide)
- Apply Fast Red solution (50µl) and observe development by light microscope.
- Wash with distilled water and **proceed rapidly** with counterstain (Harris hematoxylin at 50% dilution in distilled water for 3 dips), dehydration and coverslipping.

Emission date: 23/08/2016

Manufacturer: **CENBIMO S.L.**
C/Doctor Iglesias Otero s-n
27004-Lugo, Spain