HISTOSONDA EBER

Referencia: CBM-0001-R5

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- 1) Gómez-Román J.J, Martínez M.N, Fernández S.L, Val-Bernal J.F. Epstein-Barr virus-associated adenocarcinomas and squamous-cell lung carcinomas. *Modern Pathology* (2009) 22, 530-537
- 2) Domínguez-Duran E, et al. Laryngeal Lymphomatoid Granulomatosis in a HIV Patient. *Enferm Infecc Microbiol Clin*. 2011. doi: 10.1016/j.eimc.2011.02.017
- 3) Koletska T, et al. Synchronous Presence of Nasopharyngeal Carcinoma and Marginal Zone (MALT-Type) B-Cell Lymphoma in the Pharynx. *Pathology Research International* (2011) doi:10.4061/2011/340763
- 4) Epstein M, Achong B, Barr Y. Morphological and biological studies on a virus in cultured lymphoblast from Burkitt's lymphoma. *J.Exp.Med.* 1965; 121:761-770.
- 5) Henle G, Henle W, Clifford P, et al. Antibodies to EB virus in Burkitt's lymphoma and control groups. *J. Nat. Cancer Inst.* 1969; 43: 1147-1157.
- 6) Komano, J.S, Maruo, K, Kurozumi, T, Oda, K, Takada. 1999. Oncogenic role of Epstein-Barr virus encoded RNAs in Burkitt's lymphoma cell line Akata. *J. Virol* 73: 9827-31.
- 7) Elaine S. Jaffe, Nancy Lee Harris, Harald Stein, James W. Vardiman. 2001. Tumors of Haematopoietic and Lymphoid Tissues (World Health Organization Classification of Tumours). IARC Press. Chapters 6,7,8.

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HISTOSONDA EBER

Reference: CBM-0001-R5

Assay of one vial for 5 individual reactions. 65µl for reaction.

Product classification:
*EU countries
For In Vitro Diagnostics

INTRODUCTION

The Epstein-Barr virus is a member of the gamma-herpes viruses (HHV-4). It is a linear 184,000 base pairs double stranded DNA virus and was the first virus to be discovered as oncogenic (1).

Primoinfection by this virus can show signs of a slight viral infection or it can present as Infectious Mononucleosis. The most common target cells for the Epstein-Barr virus are the B lymphocyte and nasopharyngeal epithelial cell. The Epstein-Barr virus massively infects the human population and sero-epidemiological studies show that 90% of adults have been infected by this virus (2).

Latently infected B lymphocytes express abundantly (10^4 - 10^5 copies), among other genes, a short non-polyadenylated chain of RNA that does not transduce to a protein, consisting of two fragments known as EBER 1 and EBER 2. The expression of EBER (Epstein-Barr virus encoded RNAs) is nuclear. Although the function of EBER is unknown, it is believed that it may play a role in virus-produced oncogenesis (3).

There are numerous human tumors associated with EBV, ranging from non-differentiated nasopharyngeal carcinoma to African Burkitt's lymphoma, Hodgkin's disease Mixed cellularity, some B, T and NK lymphomas, as well as in lymphoproliferative processes associated with immunodeficiency(4).

INTENDED USE

For use in In Vitro Diagnosis. Histosonda EBER is intended for the detection of EBER 1+2 RNA in cells infected by the Epstein-Barr virus, both in reactive and tumoral cells. Its use is intended for the aforementioned tumors, as well as in lymph nodes of patients with atypical reactive pathology as produced by Infectious Mononucleosis. This detection is carried out by in situ hybridization in paraffin-embedded, formalin-fixed histological sections.

WARNINGS AND PRECAUTIONS

Histosonda EBER has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel.

In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

This product contains sodium azide as a preservative (NaN_3). Sodium azide is below the threshold level of 0,1% (w/v) of the total, and therefore it is not considered to be a hazardous substance. Small quantities may be eliminated in the waste water system.

COMPONENTS

1 tube with the necessary amount of Histosonda EBER for developing 5 single reactions.

Histosonda EBER consists of a fragment of single stranded DNA of 527 nucleotides, complementary to expressed EBER. This probe detects the EBER gene in full, including EBER 1 and EBER 2. The DNA of the probe contains nucleotides that have been labeled with Digoxigenin.

STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

INTERPRETATION OF RESULTS

Samples in which EBER expression is observed will show a brownish color in cell nucleus, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

ASSAY LIMITATIONS

Histosonda EBER has been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results. The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

BIBLIOGRAPHY

- 1) Gómez-Román J.J, Martínez M.N, Fernández S.L, Val-Bernal J.F. Epstein-Barr virus-associated adenocarcinomas and squamous-cell lung carcinomas. *Modern Pathology* (2009) 22, 530-537
- 2) Domínguez-Duran E, et al. Laryngeal Lymphomatoid Granulomatosis in a HIV Patient. *Enferm Infecc Microbiol Clin*. 2011. doi: 10.1016/j.eimc.2011.02.017
- 3) Koletska T, et al. Synchronous Presence of Nasopharyngeal Carcinoma and Marginal Zone (MALT-Type) B-Cell Lymphoma in the Pharynx. *Pathology Research International* (2011) doi:10.4061/2011/340763
- 4) Epstein M, Achong B, Barr Y. Morphological and biological studies on a virus in cultured lymphoblast from Burkitt's lymphoma. *J.Exp.Med.* 1965; 121:761-770.
- 5) Henle G, Henle W, Clifford P, et al. Antibodies to EB virus in Burkitt's lymphoma and control groups. *J. Nat. Cancer Inst.* 1969; 43: 1147-1157.
- 6) Komano J.S. Maruo K. Kurozumi T. Oda K. Takada. 1999. Oncogenic role of Epstein-Barr virus encoded RNAs in Burkitt's lymphoma cell line Akata. *J. Virol* 73: 9827-31.
- 7) Elaine S. Jaffe, Nancy Lee Harris, Harald Stein, James W. Vardiman. 2001. Tumors of Haematopoietic and Lymphoid Tissues (World Health Organization Classification of Tumours). IARC Press. Chapters 6,7,8.

PROCEDURE

For Digoxigenin- labeled probes

BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

HistoSonda Protocol

Needed accessories for this protocol and not provided with the product
Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete Histosonda Technical Manual at the email address: techserv@cenbimo.com

Before starting pre-heat the humid chamber to 62°C.

1a. Deparaffinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
 - a. Xylene: 10 mins
 - b. Xylene: 5mins
 - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
 - d. Alcohol 96%: 1min X 3
 - e. Methanol containing 0.3% H_2O_2 or 3% H_2O_2 only: 5mins
3. Wash well with distilled water, leave standing in water.

1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

1. Place slides in boiling distilled water for **30 seconds** and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033 mg/ml in PBS pH 7.4.
2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
4. Wash well with water.
5. Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment!**

3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.

4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.!

1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS pH 7.4 agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS, agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) following the manufacturer's instructions.
8. Wash with distilled water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in distilled water.
10. Wash with distilled water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: CENBIMO S.L.
C/Doctor Iglesias Otero s/n
27004-Lugo, Spain