## HISTOSONDA GHRL (Ghrelin/Obestatin)

Referencia: CBM-0021-R5

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HISTOSONDA GHRL (Ghrelin/Obestatin)

Reference: CBM-0021-R5

Assay for 5 individual reactions in one vial. 65µl for reaction.

Product classification:

\*EU countries  
\*For In Vitro Diagnostics\*

## INTRODUCTION

Ghrelin is an endogenous ligand for the secretagogue receptor of the pituitary growth hormone. Neuroendocrine tumors producing Ghrelin have been described as clinically developing with acromegaly.

The notable functions of Ghrelin, amongst others are its orexigenic action (stimulation of the ingestion of alimentation) and its stimulus of hydrochloric acid production by the stomach parietal cells.

Many studies have related some obese patients with a seric increase of Ghrelin. Some publications have demonstrated that Ghrelin plays a part in inhibiting the production of leptin.

Ghrelin is produced by neuroendocrine cells that are preferably situated in the hydrochloric acid producing glands of the gastric body.

Some authors have also described cells that produce Ghrelin in other anatomic localizations although an exhaustive "in situ" hybridization map has not yet been made.

Cenbimo has found cases of hyperplasia and tumors of gastric neuroendocrine cells that produce Ghrelin (data not yet published).

Histosonda GHRL consists of a fragment of single stranded DNA, 241 nucleotides in length targeted against Ghrelin RNA.

The use of this probe is the detection of cells producing this hormone and their tumors.

## INTENDED USE

For use in In Vitro Diagnosis.

Histosonda GHRL is useful for the detection of cells producing this hormone and their tumors.

## WARNINGS AND PRECAUTIONS

Histosonda GHRL has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel. In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

## COMPONENTS

1 tube with the necessary amount of Histosonda GHRL for developing 5 single reactions  
Histosonda GHRL consists of a segment of single-stranded DNA complementary to expressed RNA with a length of 241 nucleotides. This probe has been labeled with Digoxigenin.

## STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

## SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

## INTERPRETATION OF RESULTS

Samples in which Ghrelin expression is observed will show a brownish color in the cell cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

## ASSAY LIMITATIONS

Histosonda GHRL has been optimized to detect Ghrelin RNA expression in formalin-fixed, paraffin-

embedded tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

## PROCEDURE

For Digoxigenin- labeled probes

### BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

### HistoSonda Protocol

#### Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete **HistoSonda Technical Manual at the email address: techserv@cenbimo.com**

Before starting pre-heat the humid chamber to 62°C.

#### 1a. Deparaffinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
  - a. Xylene: 10 mins
  - b. Xylene: 5mins
  - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
  - d. Alcohol 96%: 1min X 3
  - e. Methanol containing 0.3% H<sub>2</sub>O<sub>2</sub> or 3% H<sub>2</sub>O<sub>2</sub> only: 5mins
3. Wash well with distilled water, leave standing in water.

#### 1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally **bone marrow and gastric tissues**).

1. Place slides in boiling distilled water for **30 seconds** and immediately transfer to distilled water at room temp.

! **Important:** excess heat treatment will result in background staining. Remove slides immediately. !

#### 2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033 mg/ml in PBS pH 7.4.
2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
4. Wash well with distilled water.
5. Transfer to PBS pH 7.4 for 2mins.

! **(Important:** if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment** !

#### 3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is com-

pletely covered and avoiding air bubbles.

3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.

#### 4. Washing the probe

! **The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.**!

1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

#### 5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS pH 7.4 agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS pH 7.4, agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
8. Wash with distilled water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in distilled water.
10. Wash with distilled water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: **CENBIMO S.L.**  
C/Doctor Iglesias Otero s-n  
27004-Lugo, Spain