## HISTOSONDA GLUCAGÓN

Referencia: CBM-0012-R5

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## HISTOSONDA GLUCAGON

Reference: CBM-0012-R5

Assay for 5 individual reactions in one vial. 65µl for reaction.

Product classification:  
\*EU countries  
"For In Vitro Diagnostics"

## INTRODUCTION

Insulin is produced by the pancreatic  $\beta$  cells in the islets of Langerhans. Between approximately 2/3 and 3/4 of adult pancreatic cells are producers of insulin. In the normal adult pancreas, the islets of Langerhans do not exceed a diameter of 300µm. In inflammatory or fibrotic processes the islets can be enlarged and can also be aspherical. When the islets exceed a diameter of 500µm they are referred to as "micro adenomas" and the possibility of genetic diseases such as MEN 1 (Multihormonal neoplasia 1) must be considered. In non tumoral islets the insulin producing cells are uniformly distributed over the entire islets and are very abundant. GLUCAGON producing cells have the tendency to be distributed peripherally although they can also be found in the interior of the islets. SOMATOSTATIN producing cells are the most scarce and do not demonstrate any preference of localization. Some somatostatin producing cells can have a fusiform morphology.

We rarely encounter GASTRIN and VIP (Vasoactive intestinal peptide) producing cells in the pancreatic islets and for this reason these hormones are considered ectopic to the pancreas. Although some patients have pancreatic gastrinomas and vipomas, their localization most frequently is in the duodenum.

It is more probable that benign pancreatic endocrine tumors clinically present a functional syndrome due to the peptide that they produce while malignant tumors will only rarely display a functional syndrome.

Both benign and malignant pancreatic endocrine tumors show a diversity of cells producing diverse peptides including insulin, glucagon and somatostatin, but they frequently express one of these peptides in a greater quantity than normal. Non tumoral islets of Langerhans will express this peptide at a much lesser intensity probably due to the inhibitory feedback phenomenon.

In some cases of pancreatic non-functional endocrine tumors, these tumors have been immunohistochemically found to have an absence of a protein although the RNA for this protein is present. Therefore, in situ hybridization is a much more useful technique for studying this kind of tumor and their metastasis.

## INTENDED USE

For use in In Vitro Diagnosis.

Histosonda Glucagon is useful for the study and classification of neuroendocrine tumors of the pancreas and digestive apparatus and their metastasis even in an absence of a clinically evident Endocrine Syndrome.

## WARNINGS AND PRECAUTIONS

Histosonda Glucagon has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel.

In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

This product contains sodium azide as a preservative (NaN<sub>3</sub>). Sodium azide is below the threshold level of 0,1% (w/v) of the total, and therefore it is not considered to be a hazardous substance. Small quantities may be eliminated in the waste water system.

## COMPONENTS

1 tube with the necessary amount of Histosonda Glucagon for developing 5 single reactions.

Histosonda Glucagon consists of a fragment of single-stranded DNA with a length of 414 nucleotides, complementary to expressed RNA. This probe has been labeled with Digoxigenin.

## STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

## SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

## INTERPRETATION OF RESULTS

Samples in which glucagon expression is observed will show a brownish color in the cell cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

## ASSAY LIMITATIONS

Histosonda Glucagon has been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

## PROCEDURE

For Digoxigenin-labeled probes

### BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

### HistoSonda Protocol

#### Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete Histosonda Technical Manual at the email address: [techserv@cenbimo.com](mailto:techserv@cenbimo.com)

Before starting pre-heat the humid chamber to 62°C.

#### 1a. Deparaffinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
  - a. Xylene: 10 mins
  - b. Xylene: 5mins
  - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
  - d. Alcohol 96%: 1min X 3
  - e. Methanol containing 0.3% H<sub>2</sub>O<sub>2</sub> or 3% H<sub>2</sub>O<sub>2</sub> only: 5mins
3. Wash well with distilled water, leave standing in water.

#### 1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

1. Place slides in boiling distilled water for 30 seconds and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

#### 2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033 mg/ml in PBS pH 7.4.
2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
4. Wash well with distilled water.

5. Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment) !

### 3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.
5. If no suitable humid chamber is available drop the 65µl of probe over a 24 X 50mm coverslip and place over the tissue section. Seal the edges of the coverslip with rubber cement. Incubate in a 62°C incubator or over a hot plate for 1 hour. After incubation carefully remove the rubber cement and coverslip and continue the protocol.

### 4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.!

1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

### 5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS, pH 7.4 agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS, agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
8. Wash with distilled water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in distilled water.
10. Wash with distilled water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: CENBIMO S.L.  
C/Doctor Iglesias Otero s-n  
27004-Lugo, Spain