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Assay of 4 vials (Kappa Light Chain or Lambda Light Chain) for 20 individual reactions. 65µl for reactions

Product classification:
*EU countries
"For In Vitro Diagnostics"

INTRODUCTION

Immunologically competent B lymphocytes develop from a lymphoid stem cell until acquiring a receptor consisting of a single immunoglobulin molecule specific for a determined antigen. These immunoglobulin molecules are formed by two identical heavy chains and by two identical light chains (kappa or lambda). In humans, in normal or reactive lymphoid tissues, approximately 60% of the B lymphocytes bear kappa light chains and the remaining 40%, lambda light chains. In B lymphocyte tumors, there is generally only one clone that proliferates, be it kappa or lambda, in a monoclonal way. In humans the gene that codes for kappa light chains is located on chromosome 2, and the gene that codes for lambda light chains is located on chromosome 22. The variable regions for genetic recombination are located in the first 5' half of these genes, while the constant regions are located in the 3' half. There is only one functional gene for the constant region of kappa light region, and four functional genes for the constant region of lambda light chains (CA1, CA2, CA3, and CA7). Detection of monoclonality is one of the most important tools for differentiating B lymphoid tumors from reactive processes. In situ hybridization technique offers an important advantage over immunohistochemistry, as it virtually lacks background, and allows a clean and sharp viewing of the histological preparation. It is also useful to differentiate cells that have absorbed immunoglobulins, and are therefore detectable by immunohistochemistry, but in fact do not produce immunoglobulin, as occurs with the Reed-Sternberg cells of Hodgkin's disease.

INTENDED USE

For use in In Vitro Diagnosis.
The Histosondas of the Immunoglobulin Light Chains are useful for the study of monoclonality in lymphoid tumors, lymphoproliferative syndromes, myelomas and for the study of immunodeficiency associated lymphoproliferative syndromes.

WARNINGS AND PRECAUTIONS

The Histosondas of the Immunoglobulin Light Chains have been designed for their professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel.
In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.
This product contains sodium azide as a preservative (NaN₃). Sodium azide is below the threshold level of 0.1% (w/v) of the total, and therefore it is not considered to be a hazardous substance. Small quantities may be eliminated in the waste water system.

COMPONENTS

4 tubes with the necessary amount of Histosonda Kappa or Lambda (as indicated on the label) for developing 20 single reactions.
The Histosondas of the Immunoglobulin Light Chains consist of a fragment of single-stranded DNA with a length of between 150 and 186 nucleotides, complementary to expressed RNA. These probes have been labeled with Digoxigenin.

STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

INTERPRETATION OF RESULTS

Samples in which immunoglobulin light chains expression is observed will show a brownish color in the cell

cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls for expression of this gene.

ASSAY LIMITATIONS

The Histosondas of the Immunoglobulin Light Chains have been optimized to detect RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques.
The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results. The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.
In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

BIBLIOGRAPHY

1. Peter J. Delves And Ivan M. Roitt: Immunoglobulin genes. In ENCYCLOPEDIA OF IMMUNOLOGY. Page 1323. Second Edition, ACADEMIC PRESS LIMITED (1998)
2. Weiss LM, Movahed LA, Chen YY, Shin SS, Stroup RM, Bui N, Estress P, Bindl JM. Detection of immunoglobulin light-chain mRNA in lymphoid tissues using a practical in situ hybridization method. Am J Pathol (1990) Oct; 137(4):979-88.
3. Ruprai AK, Pringle JH, Angel CA, Kind CN, Lauder I. Localization of immunoglobulin light chain mRNA expression in Hodgkin's disease by in situ hybridization. J Pathol.1991 May; 164(1): 37-40

PROCEDURE

For Digoxigenin- labeled probes

BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.
The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

HistoSonda Protocol

Needed accessories for this protocol and not provided with the product
Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete HistoSonda Technical Manual at the email address: techserv@cenbimo.com

Before starting pre-heat the humid chamber to 62°C.

1a. Deparaffinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
 - a. Xylene: 10 mins
 - b. Xylene: 5mins
 - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
 - d. Alcohol 96%: 1min X 3
 - e. Methanol containing 0.3% H₂O₂ or 3% H₂O₂ only: 5mins
3. Wash well with distilled water, leave standing in water.

1b. Inhibition of unspecific DNA binding (optional)
Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

1. Place slides in boiling distilled water for **30 seconds** and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033 mg/ml in PBS pH 7.4.

2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temp.
4. Wash well with distilled water.
5. Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment**) !

3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.

4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.!

1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS pH 7.4 agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS, agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
8. Wash with water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in distilled water.
10. Wash with distilled water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: CENBIMO S.L.
C/Doctor Iglesias Otero s-n
27004-Lugo, Spain