HISTOSONDA LISOZIMA

Referencia: CBM-0022-R5

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HISTOSONDA LYSOZYME

Reference: CBM-0022-R5

Assay of one vial for 5 individual reactions. 65µl for reaction.

Product classification:

*EU countries
For In Vitro Diagnostics

INTRODUCTION

Lysozyme, also known as Muramidase, catalyses the hydrolysis of some mucopolysaccharides in bacterial cell walls, specifically hydrolyzing between acetyl/muramic acid and acetylglucosamine.

This enzyme is found in a variety of human cells but principally in tissular histiocytes, circulating granulocytes, monocytes and in glandular saliva and lacrimal cells.

The quantity of lysozyme produced by monocytes is greatly superior to that of polymorphonuclear neutrophil granulocytes.

Various point mutations in the lysozyme gene have been associated with non-neuropathic visceral amyloidosis.

Histosonda Lysozyme consists of a fragment of single stranded DNA 358 nucleotides long that is targeted against lysozyme RNA.

The primary use of this probe is for the classification and study of histiocytic cells (especially when they are activated and functioning) and the detection of monocytes in the bone marrow.

Due to the fact that in decalcified, paraffin embedded bone marrow biopsies it is not possible to clearly see the granulocytic cytoplasmic granules it is frequently very difficult to distinguish between myelocytes and monocytes.

Histosonda Lysozyme permits a clean and clear differentiation of the bone marrow monocytic cells owing to the augmented transcription of the gene producing lysozyme. This permits the clear differentiation of chronic myelomonocytic leukemia from other types of chronic granulocytic leukemias.

INTENDED USE

For use in In Vitro Diagnosis.

Histosonda Lysozyme is useful for the classification and study of histiocytic cells (especially when they are activated and functioning) and the detection of monocytes in the bone marrow.

WARNINGS AND PRECAUTIONS

Histosonda Lysozyme has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel. In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results. This product contains sodium azide as a preservative (NaN₃). Sodium azide is below the threshold level of 0,1% (w/v) of the total, and therefore it is not considered to be a hazardous substance. Small quantities may be eliminated in the waste water system.

COMPONENTS

1 tube with the necessary amount of Histosonda Lysozyme.

Histosonda Lysozyme consists of a segment of single-stranded DNA complementary to expressed RNA with a length of 358 nucleotides. This probe has been labeled with Digoxigenin.

STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

INTERPRETATION OF RESULTS

Samples in which Lysozyme expression is observed will show a brownish color in the cell cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

ASSAY LIMITATIONS

Histosonda Lysozyme has been optimized to detect Lysozyme RNA expression in formalin-fixed, paraffin-embedded tissues. Its use is not recommended for other types of samples or preparation techniques. The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

PROCEDURE

For Digoxigenin- labeled probes

BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process. The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

HistoSonda Protocol

Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP commercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete HistoSonda Technical Manual at the email address: techserv@cenbimo.com

Before starting pre-heat the humid chamber to 62°C.

1a. Deparaffinization

- Heat slides at 62°C for 10mins using a hot plate or incubator.
- Immerse slides in:
 - Xylene: 10 mins
 - Xylene: 5mins
 - Absolute alcohol (ethanol or isopropanol): 1min X 2
 - Alcohol 96%: 1min X 3
 - Methanol containing 0.3% H₂O₂ or 3% H₂O₂ only: 5mins
- Wash well with distilled water, leave standing in water.

1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

- Place slides in boiling distilled water for 30 seconds and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

2. Deproteinization

- Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033 mg/ml in PBS pH 7.4.
- Remove excess water from individual slides with tissue paper.
- Cover tissue sections with the 200µl of Proteinase K and incubate in a humid chamber for exactly 10mins at room temperature.

- Wash well with distilled water.
- Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment!**)

3. Incubation with the probe

- Remove excess buffer from sections as described previously.
- Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
- Place the slides in a horizontal position in a humid chamber and close well.
- Incubate at 62°C for 1hr.

4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.!

- Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
- Agitate in PBS for 5mins.

5. Revealing the probe

Protocol for manual revealing of the probe.

- Remove excess buffer from sections as before.
- Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
- Wash vigorously with PBS pH 7.4 and agitate in PBS for 1min.
- Remove excess buffer from sections.
- Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
- Wash vigorously with PBS pH 7.4 and agitate in PBS 1min.
- Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
- Wash with distilled water.
- Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in distilled water.
- Wash with distilled water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: CENBIMO S.L.
C/Doctor Iglesias Otero s-n
27004-Lugo, Spain