## HISTOSONDA MMP1 (MATRIX METALOPROTEINASA 1)

Referencia: CBM-0028-R5

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## HISTOSONDA MMP1 (MATRIX METALLOPROTEINASE 1)

Reference: CBM-0028-R5

Assay of one vial for 5 individual reactions. 65µl for reaction.

### Product classification:

\*EU countries  
"For In Vitro Diagnostics"

### INTRODUCTION

The extracellular matrix is composed of structural proteins (fibronectin, laminin, fibrous collagens etc) and non-structural proteins.

The Matrix Metalloproteinases are non-structural extracellular matrix proteins that function as proteases. They have been previously known by many names such as collagenases, gelatinases, matrilysin, stromelysin etc.

The Matrix Metalloproteinases function in remodelling of the extracellular matrix and have been implicated in processes such as the tissular reorganisation in stress or injury and also in different phases of embryonic development.

Recently the Matrix Metalloproteinases have been implicated in progression, infiltration and metastasis of malignant epithelial tumours, forming part of the complex tumoural-stromal relationship.

In situ hybridisation allows RNA production of various metalloproteinases to be visualised in peritumoural stromal cells whose pattern is dependent on tumoural type and tissue (apart from MMP7 which is produced by tumoural epithelial cells).

### INTENDED USE

For use in In Vitro Diagnosis.  
Histosonda MMP1 is useful for the detection of MMP1-producing cells in tumoural and non tumoural tissues.

### WARNINGS AND PRECAUTIONS

Histosonda MMP1 has been designed for professional use in In Vitro Diagnosis and must be manipulated by qualified and accordingly trained personnel. In order to obtain the best results, the instructions contained in the manual must be followed. Any change to the indicated temperatures, times or any other step of the process can lead to poor results.

This product contains sodium azide as a preservative (NaN<sub>3</sub>). Sodium azide is below the threshold level of 0.1% (w/v) of the total, and therefore it is not considered to be a hazardous substance. Small quantities may be eliminated in the waste water system.

### COMPONENTS

1 tube with the necessary amount of Histosonda MMP1 for developing 5 single reactions.  
Histosonda MMP1 consists of two fragments of single-stranded DNA complementary to expressed RNA with lengths of 699 and 941 nucleotides. This probe has been labeled with Digoxigenin.

### STORAGE CONDITIONS

Supplied reagents are stored at 3-7°C until expiration date. Do not use after expiration date.

### SAMPLES

Block sections of formalin fixed paraffin embedded tissues. Sections of 4-6 micrometers in width are sufficient to conduct the study.

### INTERPRETATION OF RESULTS

Samples in which MMP1 expression is observed will show a brownish color in the cell cytoplasm, which will contrast over the blue-violet background given by hematoxylin staining. The pathologist will evaluate the results according to their experience, drawing conclusions from the staining of the sample, in parallel with the staining observed in the positive and negative controls.

### ASSAY LIMITATIONS

Histosonda MMP1 has been optimized to detect RNA expression in formalin-fixed, paraffin-embedded

tissues. Its use is not recommended for other types of samples or preparation techniques.

The correct operation of this product has been validated using the protocols indicated in the instructions manual. The use of other procedures or the modification of the recommended protocols may lead to erroneous results.

The results from this assay must be evaluated by the pathologist in combination with the rest of available patient clinical data.

In order to obtain optimal and reproducible results it is important to rigorously maintain the time and temperature conditions indicated in the procedure.

### BIBLIOGRAPHY

1. Egeblad M, Werb Z. New functions for the matrix metalloproteinases in cancer progression. *Nat Rev Cancer*. 2002 Mar;2(3):161-74.
2. Kitoh T, Yanai H, Saitoh Y, Nakamura Y, Matsubara Y, Kitoh H, et al. Increased expression of matrix metalloproteinase-7 in invasive early gastric cancer. *J Gastroenterol*. 2004;39(5):434-40.
3. Zeng ZS, Shu WP, Cohen AM, Guillem JG. Matrix metalloproteinase-7 expression in colorectal cancer liver metastases: evidence for involvement of MMP-7 activation in human cancer metastases. *Clin Cancer Res*. 2002 Jan;8(1):144-8.
4. Duffy MJ, Maguire TM, Hill A, McDermott E, O'Higgins N. Metalloproteinases: role in breast carcinogenesis, invasion and metastasis. *Breast Cancer Res*. 2000;2(4):252-7.
5. Overall CM, Lopez-Otin C. Strategies for MMP inhibition in cancer: innovations for the post-trial era. *Nat Rev Cancer*. 2002 Sep;2(9):657-72.

### PROCEDURE

For Digoxigenin- labeled probes

### BASIS OF THE METHOD

Chromogenic In Situ Hybridization (CISH) is a technique used to determine the presence of a DNA or RNA sequence or to study gene expression, as well as for the simultaneous evaluation of tissue morphology by white light microscopy. A labeled probe of known sequence hybridizes with its target in the tissue being studied. This hybridization is then detected by an immunohistochemical process.

The Histosondas targets are messenger RNA and they are therefore designed to detect and visualise gene expression.

### HistoSonda Protocol

#### Needed accessories for this protocol and not provided with the product

Proteinase K, Anti-Digoxin Antibody, polímero comercial anti-mouse HRP comercial polymer, Diaminobenzidine (DAB) solution.

For further important details including FAQs and Troubleshooting please ask for the complete HistoSonda Technical Manual at the email address: [techserv@cenbimo.com](mailto:techserv@cenbimo.com)

Before starting pre-heat the humid chamber to 62°C.

#### 1a. Deparaffinization

1. Heat slides at 62°C for 10mins using a hot plate or incubator.
2. Immerse slides in:
  - a. Xylene: 10 mins
  - b. Xylene: 5mins
  - c. Absolute alcohol (ethanol or isopropanol): 1min X 2
  - d. Alcohol 96%: 1min X 3
  - e. Methanol containing 0.3% H<sub>2</sub>O<sub>2</sub> or 3% H<sub>2</sub>O<sub>2</sub> only: 5mins
3. Wash well with distilled water, leave standing in water.

#### 1b. Inhibition of unspecific DNA binding (optional)

Generally this step is not necessary unless the tissue to be hybridized contains a great quantity of polymorphonuclear eosinophil leukocytes (generally bone marrow and gastric tissues).

1. Place slides in boiling distilled water for 30 seconds and immediately transfer to distilled water at room temp.

! Important: excess heat treatment will result in background staining. Remove slides immediately. !

#### 2. Deproteinization

1. Prepare a dilution of any commercial Proteinase K at a work concentration of 0,033 mg/ml in PBS pH 7.4.
2. Remove excess water from individual slides with tissue paper.
3. Cover tissue sections with the 200µl of Proteinase K and incubate in a humid

chamber for exactly 10mins at room temp.

4. Wash well with distilled water.
5. Transfer to PBS pH 7.4 for 2mins.

! (Important: if the slides have been boiled in the microwave previously only use Proteinase K for 5 minutes as the heat treatment will have significantly sensitized the tissue sections. **Bone marrows require a Proteinase K digestion of 20min at 55°C after heat treatment!** !

#### 3. Incubation with the probe

1. Remove excess buffer from sections as described previously.
2. Add 65µl of probe to the tissue section ensuring the entire section is completely covered and avoiding air bubbles.
3. Place the slides in a horizontal position in a humid chamber and close well.
4. Incubate at 62°C for 1hr.

#### 4. Washing the probe

! The probe solution is very viscous. It is very important to wash the slides well to ensure the probe solution is completely removed.!


1. Wash the surface of the section vigorously with PBS pH 7.4 to remove the probe.
2. Agitate in PBS for 5mins.

#### 5. Revealing the probe

Protocol for manual revealing of the probe.

1. Remove excess buffer from sections as before.
2. Cover sections with 100µl of Anti-Digoxin (not provided with this product) following manufacturer's instructions and incubate in a humid chamber for 30mins at room temp.
3. Wash vigorously with PBS pH 7.4 and agitate in PBS for 1min.
4. Remove excess buffer from sections.
5. Drop commercial HRP polymer (not provided with this product) that binds the primary antibody used before over the sections (enough to cover the tissue 50-100µl) and incubate in a humid chamber following the manufacturer's instructions.
6. Wash vigorously with PBS pH 7.4 and agitate in PBS 1min.
7. Remove excess buffer from sections and apply commercial Diaminobenzidine (DAB) (not provided with this product) following the manufacturer's instructions.
8. Wash with distilled water.
9. Counter stain the sections briefly (2-3secs) with Harris hematoxylin diluted 50% in distilled water.
10. Wash with distilled water, dehydrate and cover slip following normal laboratory protocols.

Emission date: 23/08/2016

Manufacturer: CENBIMO S.L.  
C/Doctor Iglesias Otero s-n  
27004-Lugo, Spain